

## OOSAFE® DISINFECTANT FOR CO2 INCUBATORS AND LAMINAR FLOW HOODS

Safe disinfection for CO<sub>2</sub> incubators, laminar flow hoods, glass, plastic, metal, Hi-Macs surfaces, equipments, ultrasound probes and hospital furnitures.



- Scientifically Proven Safety: I-Human sperm survival bioassay to examine toxicity of a new clinical laboratory equipment disinfectant (Poster at ASRM, 2010).
- 2-Changing old trends in IVF laboratory cleaning and disinfection (Study from MONASH University Australia, 2012).
- MEA Tested: Each lot is MEA tested for IVF laboratory use.
- Long Term Disinfection: Unlike volatile disinfectants (which evaporates after application), it keeps applied areas disinfected for a long time.
- No VOC: No alcohol, no smell, no VOC release.
- Long Shelf Life: 3 years after manufacturing.
- No Waste: Can be used till last drop until expiry date.
- Easy to Use: Ready to use solution, no washing or rinsing after application.
- Fast Disinfection: Complete disinfection after I5 minutes of application.
- Quality: Each bottle has security seal which assures safe transportation and quality controlled product. I liter disinfectant is packed in individual box together with spray.
- Microbiological Efficacy: Effective against hepatitis B, HIV, rota virus within I minute, mycobacterium within 5 minutes, bacteria, fungi (Candida), influenza A virus (H5NI/HINI) within I5 minutes.
- •Usage: Can be used even in laboratories which do not have enough ventilation, no VOC release after application.
- Application: Apply evenly on the surface using spray or wipe, wait I5 minutes
  to dry out for the best result. It is recommended to wipe the surface with a
  humid sterile sponge or cloth after I5 minutes to clean surface against dead
  organisms.
- Active Ingredients: Quaternary ammonium compound: Benzyl-alkyldimethyl chloride I.3 mg/g.
- OODW-70

OOSAFE® DISINFECTION WIPES

- Oosafe<sup>®</sup> Disinfection Wipes, 70 wipes/bucket
- High Quality Wipes: 260x300 mm, 40 g/m³, 100 % PP



OODIH-0050

 $\begin{array}{c} {\sf Oosafe}^{\rm e} \, {\sf Disinfectant} \, {\sf for} \, {\sf CO}_2 \\ {\sf Incubators} \, {\sf and} \, {\sf Laminar} \, {\sf Flow} \, {\sf Hoods} \\ {\sf 50} \, \, {\sf ml} \, \, {\sf with} \, {\sf spray} \end{array}$ 



OODIH-I000

Oosafe® Disinfectant for CO<sub>2</sub> Incubators and Laminar Flow Hoods I liter with spray



OODIH-5000

Oosafe® Disinfectant for CO<sub>2</sub> Incubators and Laminar Flow Hoods 5 liters refill